

Portable Ultrasonic Cleaner

NXPC-B Series



NXPC-B5020S



NXPC-B5020B

Model	NXPC-B2005S(B)	NXPC-B3010S(B)	NXPC-S5020S(B)
Capacity	5.7 Liter	10 Liter	20 Liter
Tank Size (WxDxH)	295x160x160	295x240x160	500x290x160
Overall Size (WxDxH)	325x230x320	325x330x320	530x330x297
Temp. Range	0 ~70 °C		
Timer Range	1 ~ 99 min, 4 hour		
Transducer	BLT Type		
Ultrasonic Frequency	40 KHz		
Input Power	AC220V 50/60 Hz		
Body Material	SUS-304(S) or SUS-430 Planting(B)		
Control Method	Microcontroller Control Method		
Switching Method	FET Switching		
Ultrasonic Power (Max)	200W	300W	500W
Heater Power	150W	190W	320W
Weight	9 kg	11 kg	15 kg