



M-LAB

FURNITURE LABORATORY



M-LAB





MAJOR PRODUCTS

For Basic Equipment in Lab



WWW.MAJOR-LAB.COM

M-LAB DESIGN



CONTRACTOR



MAJOR PRODUCTS

LABORATORY BENCH SET

Laboratory Bench Set is the standard laboratory desk from factory can set up in area. You can rebuild or change the cupboard or move to new location a whole or any part (MOBILE SYSTEM)




SIRENA WORKING AREA

Molded from special lab grade polyester resins. Sirena product are homogeneous in nature. There are no laminates to bubble or crack and no porosity to allow liquid penetration.


EPOXY RESIN WORKING AREA

Work top shall be molded from special grade epoxy resin and non-laminate. Chemical resistant and fire retardant. Density up to 1.90 g/cm³



MAXEL PEGBOARD

Pegboard is molded from plastic, polypropylene light weight. Pegboard with heater will dry up to 12 pieces at one time.



POSTFORM WORKTOP

PHENOLIC LAMINATE IS LAMINATE WITH HIGH PRESSURE (HPL) ON TOP AND DOOR PANEL CAN BE SELECTED IN MANY KIND OF COLOR



SINK CUPBOARD

Floor of bathtub cupboard is sloping for move the dregs-catch cup. Back wall made to be sheet, Service Panel, for move to maintenance or add other system later.



GAS SYSTEM

CONTROL WITH REGULATOR AND SAFETY OPTIONS



LABORATORY OUTLET FOR HOT OR COLD WATER

- Material: High grade brass with Epoxy powder coating surface at high temperature (300 °C) and polypropylene handle.
- Visual color : According to DIN 12920 standard
- Chemical resistance : Excellent for acid base and solvent vapor
- Flow rate: 13.25 litre/min. at 30 psi
- Maximum pressure: 120 psi.
- Line: Diameter 1/2 inch. Inlet, Diameter 3/8- 1/2 inch. Outlet .




LABORATORY DRAINAGE FITTING


- Bathtub navel is 1.5-inch diameter polypropylene with cork.
- Bathtub has a 1500 ml. Hold capacity polypropylene waste tab direct jointed with bathtub navel before jointed with building waste water system.



WORKING WOOD




M-LAB




MAJOR PRODUCTS

PARTICLE BOARD




Design



EDGE BANDING



HEAT RESISTANCE



HPL LAMINATE EDGING



MELAMINE




CURVE EDGING

MOISTURE RESISTANT
CHIPBOARD WORKTOP



KNOCK DOWN WITH
QUIPKFIT



HIGH GRADE FITTING



ABS. ADJUSTABLE LEG



ADJUSTABLE LEG

MAJOR PRODUCTS