

# Portable JAC Series

## JAC -Series



### Features

Microcomputer control system enables easy control of timer, output, and temperature.

The use of Industrial BLT transducer increases its durability and power.

Ability to control the ultrasonic output intensity – High, Moderate, Low

Ability to change the ultrasonic operation mode – Normal, Pulse, Sweep

DeGAS function is added – dissolved gas elimination to improve cleaning quality.

Inflammable ABS coated surface creates clean and elegant appearance.

### Specification

Name	JAC Series Portable Ultrasonic Cleaner	
Model	JAC-3010	JAC-5020
Capacity	10 Liter	20 Liter
Tank Size (WxDxH)	295x240x160	500x290x160
Overall Size (WxDxH)	490x368x297	721x451x297
Transducer	BLT Type	
Ultrasonic Power(Max)	300W	500W
Heater Power	190W	320W
Frequency	40 KHz	
Input Power	Single-Phase AC220V 50/60Hz	
Temp. Range	Room temperature ~ 70 °C	
Time Range	1 ~ 99 min	
Body Material	ABS	
Weight	11 kg	17 kg



### Warnings

**Please carefully read for proper and safe operation of the product.**

1. Avoid cleaning solvent that contains Benzene or Alcohol to prevent inflammable ABS coated surface from melting down.
2. Please fill water up to the indicated line before you turn on the power.
3. Heater should be on only when the machine is in use.
4. Product warranty period is one year from date of purchase.

